

Phonak Paradise 90 RT

Tests conducted with quality wireless headphones in quiet room.

Volume settings identical for all tests and measured against a calibration file.

The test files are based on the ISO 389-7:2005 international standard.

Values are averages of 3 tests each with and without HA.

N = without hearing aids

Y = with hearing aids.

	250		500		1000		2000		4000		8000	
	N	Y	N	Y	N	Y	N	Y	N	Y	N	Y
0												
10	O	O										
15	X	X										
20						O						
30			X O	X O	X O	X		O				
40							O					
50								X	O	O	O	
60							X					O
70												
80									X	X	X	X

250 and 500 – no difference in either ear.

1000 – no difference in left ear, 10 dB improvement in right.

2000 – 10 dB improvement in both ears

4000 – no difference in either ear

8000 – 10 db decline in right ear, no change in left.

Left ear at 8K cannot be measured with this test. Original audiogram measured left at 105 dB.