


Client information


Hearing Care Professional	sergo	Name	Sergey Koval
Session date	05.04.18	Date of birth	29.10.67
Client number	G00002	Age	18 - 59 years

Audiogram

RIGHT EAR



LEFT EAR



Hearing instrument settings

(R) Personalization

Experience None

(L) Personalization

Experience None

(R) Opn 2 miniRITE (105)

Serial number Not present  
Battery size 312

(L) Opn 2 miniRITE (105)

Serial number Not present  
Battery size 312

(R) Acoustics

Earpiece Power Mold  
Vent 0.8 mm

(L) Acoustics

Earpiece Power Mold  
Vent 0.8 mm

Hearing Care Professional	sergo	Client number	G00002
Session date	05.04.18	Name	Sergey Koval

**(R) Buttons**

Volume control	On
Program switch	On
Mute	On

**(L) Buttons**

Volume control	On
Program switch	On
Mute	On

**(R) Instrument setup**

Adaptation step	1
Feedback margin	Predicted
Binaural broadband	On

**(L) Instrument setup**

Adaptation step	1
Feedback margin	Predicted
Binaural broadband	On

**(R) Data logging**

Logging	On
---------	----

**(L) Data logging**

Logging	On
---------	----

**(R) Autophone**

Autophone	Off
-----------	-----

**(L) Autophone**

Autophone	Off
-----------	-----

**(R) P1: GENERAL, VAC+**

	125	250	500	625	750	1k	1.25k	1.5k	1.75k	2k	3k	4k	5k	6k
MPO	109	118	112	111	110	111	115	117	120	125	124	119	126	111
Loud	1	3	4	5	6	7	8	13	15	18	19	18	17	18
Moderate	1	6	8	9	10	13	17	22	26	28	30	29	28	27
Soft	4	11	15	17	19	23	28	34	35	40	30	29	29	27
CR	1,4	1,4	1,5	1,5	1,4	1,8	2,5	2,5	3,7	3,0	3,7	3,7	3,7	2,5
CR	1,4	1,4	1,4	1,5	1,5	1,4	1,7	1,9	1,5	1,7	1,0	1,0	1,1	1,0



Hearing Care Professional	<b>sergo</b>	Client number	<b>G00002</b>
Session date	<b>05.04.18</b>	Name	<b>Sergey Koval</b>

**(L) P1: GENERAL, VAC+**

	125	250	500	625	750	1k	1.25k	1.5k	1.75k	2k	3k	4k	5k	6k
MPO	109	118	112	111	110	111	115	117	120	125	124	119	126	111
Loud	1	3	4	5	6	7	8	13	15	18	19	18	17	18
Moderate	1	6	8	9	10	13	17	22	26	28	30	29	28	27
Soft	4	11	15	17	19	23	28	34	35	40	30	29	29	27
CR	1,4	1,4	1,5	1,5	1,4	1,8	2,5	2,5	3,7	3,0	3,7	3,7	3,7	2,5
CR	1,4	1,4	1,4	1,5	1,5	1,4	1,7	1,9	1,5	1,7	1,0	1,0	1,1	1,0

**(R) P1: GENERAL, VAC+**

**(L) P1: GENERAL, VAC+**



Hearing Care Professional	sergo	Client number	G00002
Session date	05.04.18	Name	Sergey Koval

**(R) P1: GENERAL, VAC+**

Anti-feedback system	On
Reserve gain	10 dB
Transient Noise Management	Medium
Spatial Noise Management	On
Noise reduction	On
Noise reduction - Simple	-1.5 dB
Noise reduction - Complex	-5 dB
Directionality	Open Automatic
OpenSound transition	Medium
Tinnitus sound	Off
Speech Rescue	Off

**(L) P1: GENERAL, VAC+**

Anti-feedback system	On
Reserve gain	10 dB
Transient Noise Management	Medium
Spatial Noise Management	On
Noise reduction	On
Noise reduction - Simple	-1.5 dB
Noise reduction - Complex	-5 dB
Directionality	Open Automatic
OpenSound transition	Medium
Tinnitus sound	Off
Speech Rescue	Off